

**1. Part Number:** A643-500**2. Photograph:****3. Product Information:**

Description:	Industrial Scoop
Nominal Scoop Volume:	500ml
Material of Construction:	Stainless Steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	290mm
Body Length:	185mm
Max. Body Width:	105mm
Max. Body Height:	65mm
Nominal Weight:	235g
Recommended Storage Conditions:	Dry and ambient temperature

**Note:** Due to the method of manufacture these scoops may not be suitable for hygienic applications. For hygienic applications we recommend our range of PharmaScoops®.

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Apr-24

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.