

**1. Part Number:** A674-215

**2. Photograph:**



**3. Product Information:**

Description	Instrument Tray 215 x 160mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Overall Length (including rim):	215mm
Overall Width (including rim):	160mm
Height:	18mm
Nominal Weight:	220g
Storage Conditions:	Dry, ambient conditions

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 13-Oct-20