

1. Part Number: A674-215

2. Photograph:



3. Product Information:

Description	Instrument Tray 215 x 160mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Overall Length (including rim):	215mm
Overall Width (including rim):	160mm
Height:	18mm
Nominal Weight:	220g
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25