

**1. Part Number:** A674-295

**2. Photograph:**



**3. Product Information:**

Description	Instrument Tray 295 x 197mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Overall Length (including rim):	295mm
Overall Width (including rim):	197mm
Height:	18mm
Nominal Weight:	415g
Storage Conditions:	Dry, ambient conditions

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20