

1. Part Number: A674-395

2. Photograph:



3. Product Information:

Description	Instrument Tray 395 x 275mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Overall Length (including rim):	395mm
Overall Width (including rim):	275mm
Height:	11mm
Nominal Weight:	915g
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20