

1. Part Number: A674-490

2. Photograph:



3. Product Information:

Description	Instrument Tray 490 x 320mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Overall Length (including rim):	490mm
Overall Width (including rim):	320mm
Height:	15mm
Nominal Weight:	1040g
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25