

1. Part Numbers: A684-240

2. Photograph:




3. Product Information:

Description	General Purpose Tray 240x160mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1.2L
Overall Length (including rim):	240mm
Overall Width (including rim):	160mm
Height:	50mm
Nominal Weight:	300g
Storage Conditions:	Dry, ambient conditions
Note	This product should be cleaned before first used

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25