

**1. Part Numbers:** A684-310

**2. Photograph:**



**3. Product Information:**

Description	General Purpose Tray 310x197mm
Material of Construction:	304 stainless steel
Method of Construction:	Pressed, single piece construction
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1.5L
Overall Length (including rim):	310mm
Overall Width (including rim):	197mm
Height:	47mm
Nominal Weight:	395g
Storage Conditions:	Dry, ambient conditions
Note	This product should be cleaned before first used

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25