


1. Part Number: A686-15L

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description | Lid for 15L Bucket with Spout |
| Material of Construction: | 316L stainless steel |
| Method of Construction: | Handle spot welded to top of lid Crevice free interior |
| Surface Finish: | Better than 0.5 microns Ra |
| Max. width over spout: | 380mm |
| Diameter: | 330mm |
| Nominal Weight: | 1kg |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin. |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Note: | Bucket with Spout 15L (A686-15) is sold separately to the lid. |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 8-Jan-25