


1. Part Number: A716-200

2. Photograph:



3. Product Information:

Description:	Powder Funnel 200mm
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	3.8 Litres
Top Diameter:	200mm
Overall Height:	232mm
Spout Length:	30mm
Internal Spout Diameter:	34.9mm (Spout ends with a 1.5" tri-clamp ferrule BS4825-3)
Nominal Weight:	1240g
Quantity per Pack:	1
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 8-Jan-25