


**1. Part Number:** A746-03

**2. Photograph:**



### 3. Product Information:

Description:

Jacketed Container 3L

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)

Method of Construction:

All welds ground and polished  
Crevice free interior

Surface Finish:

Better than 0.5 microns Ra

Nominal Volume:

3 Litres

Dimensions:

- A = 207mm
- B = 160mm
- C = 190mm
- D = 40mm
- E = 115mm
- F = 1.2mm
- G = 119mm

Neck / Lid Size:

4" BS4825-3

Jacket Inlet/Outlet Size:

Compatible with 1" and 1.5" tri-ferrule connections BS4825-3



BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Weight:

3.4 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions:

Dry and ambient temperature

#### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Oct-21