

1. Part Number: A764-30

2. Photograph:



3. Product Information:

Description:	Medicine Cup 30ml
Material of Construction:	304 stainless steel
Method of Construction:	Single piece construction
Surface Finish:	Better than 0.8 microns Ra
Top Diameter:	52mm
Base Diameter:	32mm
Height:	25mm
Nominal Weight:	14g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25