

1. Part Number: A776-3000

2. Photograph:



3. Product Information:

Description	GMP PharmaScoop® with Top Handle 3000ml
Nominal Volume:	3000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	420mm
Body Length:	270mm
Body Diameter:	160mm
Nominal Weight:	1.4kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 17-May-23