


1. Part Number: A776-5000

2. Photograph:



3. Product Information:

Description	GMP PharmaScoop® with Top Handle 5000ml
Nominal Volume:	5000ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	500mm
Body Length:	350mm
Body Diameter:	195mm
Nominal Weight:	1745g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25