


**1. Part Number:** A836-1000

**2. Photograph:**



**3. Product Information:**

Description:	Stainless Bottle 1000 ml
Material of Construction:	Body: 316L stainless steel Lid: 316L stainless steel Clamp: 304 stainless steel Gasket: Silicone (conforms to FDA CFR 177.2600)
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	1000ml
Brimful Volume:	1200ml
Dimensions:	A = 64 mm B = 48 mm C = 178 mm D = 102 mm
Nominal Weight (body, lid, clamp & gasket):	980 g
Recommended Storage Conditions:	Dry and ambient temperature
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin

Note: The neck of the container is fitted with a 2" diameter ferrule which conforms to BS4825-3

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 23-Feb-21