


1. Part Number: A846-1000

2. Photograph:



3. Product Information:

| | |
|---------------------------------------------|-------------------------------------------------------------------------------------------------------------------|
| Description | Stainless PharmaBottle 1000ml |
| Materials of Construction: | - Body: 316L stainless steel - Lid: 316L stainless steel - Gasket: Silicone (FDA 21 CFR 177.2600 compliant) |
| Nominal Volume: | 1000ml |
| Brimful Volume: | 1150ml |
| Method of Construction: | Crevice free interior, welds ground & polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Dimensions: | A = 69 mm B = 48 mm C = 178 mm D = 102 mm |
| Nominal Weight (body, lid, clamp & gasket): | 895g |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin. |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20