


**1. Part Number:** A846-500

**2. Photograph:**



### 3. Product Information:

Description	Stainless PharmaBottle 500ml
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Gasket: Silicone (FDA 21 CFR 177.2600 compliant)
Nominal Volume:	500ml
Brimful Volume:	650ml
Method of Construction:	Crevice free interior, welds ground & polished
Surface Finish:	Better than 0.5 microns Ra
Dimensions:	A = 69 mm B = 48 mm C = 154 mm D = 83 mm
Nominal Weight (body, lid, clamp & gasket):	690 g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20