


1. Part Number: A856-100

2. Photograph:



3. Product Information:

Description	Stainless Container with 6" Top 100L
Material of Construction:	Body: 316L stainless steel Lid: 316L stainless steel Clamp: 316L stainless steel (may be 304 stainless steel) Gasket: Silicone (conforms to FDA 21 CFR 177.2600)
Method of Construction:	Fully welded, the weld is ground and polished
Surface Finish:	Better than 0.5 microns Ra
Dimensions:	A = 167mm B = 150mm C = 670mm D = 465mm E = 500mm Wall Thickness = 2mm
Nominal Volume:	100L
Graduation Marks:	None
Pressure Rating:	None
Nominal Weight (body, lid, clamp & gasket):	20.9 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Storage Conditions:	Dry, ambient conditions

Note: The neck of the container is fitted with a 6" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25