


1. Part Number: A856-75

2. Photograph:



3. Product Information:

| | |
|---|---|
| Description | Stainless Container with 6" Top 75L |
| Material of Construction: | Body: 316L stainless steel Lid: 316L stainless steel Clamp: 316L stainless steel (may be 304 stainless steel) Gasket: Silicone (conforms to FDA 21 CFR 177.2600) |
| Method of Construction: | Fully welded, the weld is ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Dimensions: | A = 167mm B = 150mm C = 650mm D = 465mm E = 450mm Wall Thickness = 2mm |
| Nominal Volume: | 75L |
| Graduation Marks: | None |
| Pressure Rating: | None |
| Nominal Weight (body, lid, clamp & gasket): | 17.7 Kg |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Storage Conditions: | Dry, ambient conditions |

Note: The neck of the container is fitted with a 6" diameter ferrule which conforms to BS4825-3

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25