


**1. Part Number:** A886-6000

**2. Photograph:**



### 3. Product Information:

Description	Stainless Bottle GL45 - 6000ml
Materials of Construction:	- Body: 316L stainless steel - Lid: 316L stainless steel - Gasket: PTFE (FDA 21 CFR 177.1550 compliant)
Nominal Volume:	6000ml
Brimful Volume:	6070ml
Method of Construction:	Upper and lower sections welded. Interior weld left as laid but electropolished.
Surface Finish:	Electropolished (internal and external) Better than 0.5 microns Ra
Dimensions:	A = 50 mm B = 28 mm C = 290 mm D = 181 mm
Nominal Weight (body, lid & gasket):	1.48kg
BSE/TSE statement:	BSE/TSE Free (the electropolishing process removes all polishing waxes)
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 15-Dec-25