

1. Part Number: A904-165

2. Photograph:



3. Product Information:

Description	Funnel with Strainer
Material of Construction:	304 stainless steel
Side Handle:	Spot welded to cone section of funnel
Surface Finish:	Better than 0.5 microns Ra
Top Diameter:	165mm
Outer Diameter of Spout:	12mm
Overall Height:	165mm
Diameter of Holes in strainer:	4mm
Is strainer removeable?	No
Nominal Weight:	212g
Storage Conditions:	Dry, ambient conditions

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Oct-20