

1. Part Numbers: A916-1500

2. Photograph:



3. Product Information:

Description	Paddle with Holes 1500mm
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	1500mm
Diameter of Handle:	25mm
Blade Length:	300mm
Blade Width:	150mm
Nominal Weight:	2.1 Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 8-Jan-25