

1. Part Number: A976-450

2. Photograph:



3. Product Information:

Description	Angled Forceps 450mm
Material of Construction:	316L stainless steel
Method of Construction:	Spot welded and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	450mm
Nominal Weight:	274g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 16-May-25