

1. Part Number: M1000B-250

2. Photograph:



3. Product Information:

Description:

Metal Detectable Scoop 250ml - Blue

Material of Construction:

Polypropylene (PP)

- Metal detectable grade
- Only virgin material is used
- PP conforms to FDA 21 CFR 177.1520
- Blue colourant conforms to FDA 21 CFR 178.3297
- Filler conforms to FDA 21 CFR 186.1374
- Conforms to EU Regulations 10/2011
- Conforms to EC Regulations 1935/2004
- BSE/TSE free

Method of Construction:

Injection moulded as a single piece

Overall Length:

209mm

Handle Length:

104mm

Scoop Width:

81mm

Nominal Weight:

65g

Moulding Environment:

Class 100,000 Cleanroom

Packing Environment:

Class 10,000 Cleanroom

Individually Bagged?

Yes (heat sealed PE/PET bag)

Number of scoops per Pack:

4

Recommended Storage Conditions:

Dry and ambient temperature

Shelf Life:

5 years from the date of manufacture

ATEX:

This device has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Apr-19