

**1. Part Number:** M1000B-60

**2. Photograph:**



**3. Product Information:**

Description:	Metal Detectable Scoop 60ml - Blue
Material of Construction:	Polypropylene (PP) <ul style="list-style-type: none"> <li>• Metal detectable grade</li> <li>• Only virgin material is used</li> <li>• PP conforms to FDA 21 CFR 177.1520</li> <li>• Blue colourant conforms to FDA 21 CFR 178.3297</li> <li>• Filler conforms to FDA 21 CFR 186.1374</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Injection moulded as a single piece
Overall Length:	178mm
Handle Length:	106mm
Scoop Width:	49mm
Nominal Weight:	19g
Moulding Environment:	Class 100,000 Cleanroom
Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Number of scoops per Pack:	4
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare & Metal Detectable ranges please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 4-Jul-19