

1. Part Number: M1051B-200

2. Photograph:



3. Product Information:

Description:	Metal Detectable Powder Funnel 200mm - Blue
Material of Construction:	Polypropylene (PP) <ul style="list-style-type: none">• Metal detectable grade• Only virgin material is used• PP conforms to FDA 21 CFR 177.1520• Blue colourant conforms to FDA 21 CFR 178.3297• Filler conforms to FDA 21 CFR 186.1374• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Injection moulded as a single piece
Max. Diameter:	200mm
Overall Height:	165mm
OD at Bottom of Spout:	43mm
ID at Bottom of Spout:	40mm
Nominal Weight:	141g
Moulding Environment:	Class 100,000 Cleanroom

Packing Environment: Class 10,000 Cleanroom
Individually Bagged? Yes (heat sealed PE/PET bag)
Quantity per Pack: 4
Recommended Storage Conditions: Dry and ambient temperature
Shelf Life: 5 years from the date of manufacture
ATEX: This device has not been ATEX tested

4. Country of Origin: Made and packed in the UK.

5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 25-Jan-23