


 <b>SAMPLING SYSTEMS</b>	<h1>Sampler Specification</h1>	Specification Number	8053B-50
		Revision	01

Part No. 8053B-50

Description SteriWare® Ladle Blue PP

Material of Construction: Polypropylene (PP)

- FDA acceptable CFR 177.1520
- EU regulations 10/2011 compliant
- EC regulation 1935/2004 compliant



Nominal Length 260mm

Diameter of cup: 48mm

Nominal volume 50ml

Nominal Weight of Sampler: 15g


Moulding & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE bag)

Method of Sterilisation: N/A

Number of Samplers per Box: 40

Recommended Storage Conditions: Dry and ambient temperature

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 29-Sep-17
----------------------------	--	--