

	<h1>Product Specification</h1>	Specification Number	8890P-250S
		Revision	01

Part No. 8890P-250S

Description: SteriWare® Powder Bottle 250ml

Material of Construction: Polypropylene (PP) – opaque

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Nominal Volume: 250ml

Brimful Volume: 280ml

Method of Construction: Single piece, injection moulded

Height: 88mm

Lid Diameter: 76mm

Nominal Weight: 56g

Assembly & Packing Environment: Class 100,000 Cleanroom

Individually Bagged? Yes (heat sealed PE/PET bag)


Method of Sterilisation: Gamma Irradiation (25-45kGy)

Number of containers per Box: 50

BSE/TSE Free: Yes

Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade mark belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 04 April 2017
----------------------------	--	--