

	<h1>Product Specification</h1>	Specification Number	8890P-250
		Revision	01


Part No. 8890P-250  
Description: SteriWare® Powder Bottle 250ml  
Material of Construction: Polypropylene (PP) – opaque

- Conforms to FDA CFR 177.1520
- Conforms to EU Regulation 10/2011
- Conforms to EC Regulation 1935/2004



Nominal Volume: 250ml  
Brimful Volume: 280ml  
Method of Construction: Single piece, injection moulded  
Height: 88mm  
Lid Diameter: 76mm  
Nominal Weight: 56g  
Assembly & Packing Environment: Class 100,000 Cleanroom  
Individually Bagged? Yes (heat sealed PE/PET bag)  
Method of Sterilisation: N/A  
Number of containers per Box: 50  
BSE/TSE Free: Yes  
Recommended Storage Conditions: Dry and ambient temperature

SteriWare is a registered trade mark belonging to Sampling Systems Ltd.

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 04 April 2017
----------------------------	----------------------------------------------------------------------------------------------------	--------------------------------------------------------------------