


	<h1>Product Specification</h1>	Specification Number	A416-30
		Revision	01

Part No. A416-30

Description Stainless Container 30L

Materials of Construction:

- Body: 316 stainless steel
- Lid: 316 stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)



Method of Construction: All welds ground and polished
Crevice free interior

Surface Finish: Better than 0.5 microns Ra


Nominal Volume: 30 L

Dimensions

- A = 119mm
- B = 98mm
- C = 602mm
- D = 303mm

Weight: 6.90 Kg (body + lid + gasket + clamp)

Note: The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

Specification Approved By:	Signature: 	Name: Jonathan Fenner Position: Director Date: 11-Dec-17
----------------------------	--	--